

# Quick start protocol

1.

**Add** your sample in appropriate buffer to the substrate (remember controls and replicas)

2.

**Incubate** the reaction at a certain temperature with shaking

3.

**Centrifuge** (max 12,000 × g)

4.

**Transfer** the reaction product into the cuvette or other chosen format for the spectrophotometric analysis

5.

**Detect** the absorbance (404 nm for yellow, 517 nm for red, 595 nm for blue and 630 nm for green substrate)

More information can be found on our webpage: [www.glycospot.dk](http://www.glycospot.dk)

GlycoSpot 

*Enzyme screening standard*

Thorvaldsensvej 40 - 1871 Frederiksberg, Denmark